



APERTURE
Ø .39 [10]

4-40 TAPPED HOLES

▽ .09 [2.3]
4 PL


[30]
1.181

[79]
3.103

[30]
1.181

[56]
2.20


REV.	DESCRIPTION	BY	DATE
0	INITIAL CONTROL.	cmansfield	4/52016



M4x0.7 ▽ .20

[39.8]
1.565

[28.2]
1.112



DIMENSIONS: INCHES[mm]

UNSPECIFIED LIMITS DIMENSIONS ARE IN INCHES		DO NOT SCALE DRAWING		HINDS INSTRUMENTS Hinds Instruments, Inc. 7245 NE Evergreen Pkwy, Suite 150 Hillsboro, OR 97124 PHONE (503) 690-2000 FAX (503) 690-3000		
FRACTIONS	± 1/32	DATE PRINTED:		POLSnap, Overall Dimensions		
DECIMALS	.XX ± .010 .XXX ± .005	APPROVALS				DATE
HOLE DIA's	.XX ± .010 .XXX ± .003 .XXX ± .001	DRAWN				4/5/16
ANGLES	± 0° 30'	CHECKED				
BENDS	± 2°	MATERIAL		FINISH	SIZE B	
PERPENDICULARITY	.003/INCH	SAMPLES MUST BE APPROVED BY ENGINEERING DEPARTMENT PRIOR TO STARTING PRODUCTION		FINISH	FSCM NO.	
CONCENTRICITY	.003 TIR	THIS DRAWING IS THE PROPERTY OF HINDS INSTRUMENTS INC. ALL RIGHTS RESERVED.		FINISH	DWG. NO.	
				FINISH	04052016-01-OD	
				FINISH	REV. 0	
				FINISH	SCALE 1:1	
				FINISH	RELEASE DATE:	
				FINISH	SHEET 1 OF 1	